


HELENA LABORATORIES MYELOMA MODULE •WORKFLOW•

1

COLLECT PATIENT SAMPLE



**POSITIVE ID
AUDIT TRAIL**
Tube to Report



Patient info
from LIS



Automatic barcode reading with
ESH Touch or **V8 E-Class**

2

TESTING PERFORMED




ESH TOUCH
AUTO SAMPLING
AUTO DILUTIONS



SPIFE TOUCH
AUTO
PROCESSING



QS TOUCH
AUTO SCANNING
ADVANCED REPORTING



V8 E-CLASS
AUTO CE
AND SAMPLE PREP

3

PATIENT RESULTS ADDED TO DATABASE


Data from both the **QS** and the **V8** are combined in a single networked library and paired with current results for remote interpretation using the QuickScan Touch Myeloma Module.



4

PATHOLOGIST REVIEW WITH MYELOMA MODULE

**ALL HISTORICAL
DATA
DISPLAYED**



**PATHOLOGIST
REVIEWS
AND INTERPRETS**



ACCEPT

DATA QUEUED
FOR REPORTING
TO LIS

RE-EDIT

BACK TO THE TEST
BENCH FOR
ADDITIONAL TESTING



800-231-5663
www.helena.com

549-277
12/16(2)