



New! ICHOR® II Analyzer

Cat. No.	Item
5880	ICHOR II Analyzer
5882	ICHOR II E-Z Cleanser (100 mL)
5883	ICHOR II Probe Cleanser (17 mL)
5884	ICHOR II Diluent Solution (20L)
5885	ICHOR II Rinse Solution (5L)
5886	ICHOR II Lyse Reagent (500 mL)
5887	ICHOR II Tri-Level Control
5888	ICHOR II Calibrator

The ICHOR II is a lower-cost impedance cell counter for use with Plateletworks and hematology testing at the point of care. When used with Plateletworks reagents, the ICHOR II system provides a hematology profile of a standard CBC with platelet count plus platelet aggregation or inhibition in under five minutes. For ease of use, ICHOR II offers touch screen capabilities and cap piercing and supports optional use of PS/2 keyboard, barcode scanner or printer. ICHOR II utilizes reagent bottles with reagent inventories managed by the instrument.



ICHOR® Hematology Analyzer

Cat. No.	Item
ICH001	ICHOR Hematology Analyzer, 110/220V
ICH-eSAT	ICHOR eSAT Data Manager (110/220V)
IP001	ICHOR Reagent Pack
IQC-45	ICHOR MICROS 45 Hematology Control Kit
IQC-60	ICHOR MICROS 60 Hematology Control Kit
IQC-CAL	ICHOR MICROS Hematology Calibrator
ISC-MEM	ICHOR Memory SmartCard
1720	Printer Paper Pack (5 rolls)

The ICHOR hematology instrument is a compact, easy-to-use, automated impedance cell counter. The safety and ease of automated closed-tube sampling, coupled with a self-contained reagent pack and waste container, make ICHOR well-suited to point of care use.

Quality Control can be run in minutes with results stored for analysis and display.

ICHOR II ICHOR - eSAT