

Investigation of Platelet Function Disorders


Patient History & Physical



Patient Symptoms (Response)



Family History



Plateletworks[®] 
Get The Whole Platelet Picture

Helena Laboratories
Point of Care 

Plateletworks Results

Collagen	ADP	Arachidonic Acid	Interpretation
Normal	Normal	Normal	No significant effects of anti-platelet drugs or factors
Abnormal	Abnormal	Normal	Probable ADP inhibitor or anti-IIb/IIIa present
Abnormal	Abnormal	Abnormal	Probable ADP inhibitor or anti-IIb/IIIa with probable COX-1 inhibitor
Abnormal	Normal	Abnormal	Probable COX-1 inhibitor
Abnormal	Normal	Normal	Probable fibrinolysis or adhesion defect

Plateletworks Products

Cat. No.	Item
5860	Plateletworks ADP Kit (25 Tests)
5862	Plateletworks Collagen Kit (25 Tests)
5864	Plateletworks Arachidonic Acid Kit (25 Tests)
5850	Plateletworks Combo-15 Kit (15 EDTA, 15 ADP, 15 Collagen)
5852	Plateletworks Combo-25 Kit (25 ADP, 25 Collagen)
5854	Plateletworks Combo-10 Kit (10 EDTA, 10 ADP, 10 Collagen, 10 Arachidonic Acid)