

Month/Year:

Serial Number:

### V8 Nexus Preventive Maintenance Schedule

DAILY	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Empty Sample Cup Waste Drawer																															
Empty Liquid Waste Bottle																															
Fill Sample Cup Tower																															
Running Buffer Added Or Changed																															
Storage Buffer Added Or Changed																															
Maint. Buffer Added Or Changed																															
Diluent Added Or Changed																															
Clean Needle with alcohol wipe																															
Technician Initials																															
WEEKLY																															
Wipe instrument clean																															
Technician Initials																															

Reagent Changed	Lot Number	Date Changed	Tech

Reagent Changed	Lot number	Date Changed	Tech

Comments:

**Note: Duplicate this page for maintenance copies. Initial each item as required. Refer to the Operator's Manual for details. If you are uncertain about how to perform any step, please contact Helena Laboratories Electronic Customer Service Dept. (1-800-231-5663) for more information.**