

HELENA LABORATORIES

KIT SAFETY DATA SHEET

Kit Product Name: QuickGel SPE Vis Kit

Kit Catalog Number(s): 3405, 3505T

Date Prepared: 8/28/2025

Revision: 5

Kit Contents:

Component Name	Component Number
QuickGel SPE Gel	552305
Acid Blue Stain	551151
Citric Acid Destain	551959

HELENA LABORATORIES**SAFETY DATA SHEET**

DATE PREPARED: 8/28/2025

REVISION: 1

1. IDENTIFICATION**Product number:** 552575, 552305**Product identifier used on the label:** SPE Gel**Recommended use of the chemical and restrictions in use:** For In-Vitro Diagnostic use**Name, address, and telephone number of the chemical manufacturer, importer, or other responsible party:**

Helena Laboratories
1530 Lindbergh Dr.
Beaumont, TX 77707
Tel: (409) 842-3714
USA Toll Free 800-231-5663

Emergency phone number:

(409) 842-3714

2. HAZARD IDENTIFICATION**GHS Classification in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200)**

Not classified

GHS Label elements

No signal word, no hazard statement(s), no precautionary statement(s) required

Other Hazards

None known.
Handle with standard laboratory precautions.

3. COMPOSITION/INFORMATION ON INGREDIENTS**Mixtures:**

This product contains no other hazardous chemicals at or above OSHA threshold concentrations.

4. FIRST AID MEASURES

Inhalation: Remove person to fresh air. If breathing is difficult, seek immediate medical attention.

Skin contact: Wash thoroughly with soap and water. Get medical attention if irritation develops.

Eye contact: Immediately flush eyes with water for several minutes. Remove contact lenses if present. Continue rinsing.

Ingestion: If swallowed: Immediately make victim drink water. Consult a physician.

5. FIREFIGHTING MEASURES

Suitable extinguishing media: Water spray, Dry chemical, carbon dioxide, foam.

Special hazards arising from the substance or mixture: Emit toxic fumes under fire

Special protective equipment and precautions for firefighters: Wear proper protective equipment and self-contained breathing apparatus for firefighting if necessary.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment, and emergency procedures: Wear personal protective equipment. Soak up with inert absorbent material and transfer to container for proper disposal.

Methods and materials for containment and cleaning up: Do not allow material to enter drains/surface water/ground water. Ventilate the area.

7. HANDLING AND STORAGE

Wear appropriate PPE
Store in a cool, dry well-ventilated area.

Protect from physical damage.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Exposure limits (ACGIH TLV & OSHA PEL): No Data Available

Exposure controls:

Respiratory protection not require

Skin protection: Wear appropriate protective gloves and suitable protective clothing.**Eye protection:** Wear appropriate safety glasses or goggles.**Hygiene measures:** Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES**Appearance:** Solid (gel)**Odor:** None**Odor threshold:** No data available**pH:** Not applicable**Melting point/Freezing point:** No data available**Initial boiling point and range:** No data available**Flash point:** No data available**Evaporation rate:** No data available**Flammability (solid, gas):** No data available**Upper/lower flammability or explosive limits:** No data available**Vapor pressure:** No data available**Relative density:** No data available**Solubility in water:** Soluble**Partition coefficient: n-octanol/water:** No data available**Auto-ignition temperature:** No data available**Decomposition temperature:** No data available**Viscosity:** No data available

10. STABILITY AND REACTIVITY

Reactivity: Not normally reactive

Chemical stability: Stable under normal condition

Possibility of hazard reaction: Hazardous polymerization will not occur

Conditions to avoid: Avoid contact with incompatible materials and heat

Incompatible materials: Oxidizing agents

Hazard decomposition products: No Data Available

11. TOXICOLOGICAL INFORMATION

Product information: No acute toxicity information is available for this product

Potential health effects:

Skin: May cause skin irritation

Eye: May cause eye irritation

Inhalation: May be harmful if inhaled

Ingestion: May be Harmful if swallowed

Carcinogenicity: No data available

Other important toxicological hazards: To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

12. ECOLOGICAL INFORMATION

Ecotoxicity: Not expected to be hazardous to the environment

Persistence and degradability: No data available

Bio-accumulative potential: No data available

Mobility in soil: No data available

Other adverse effects: No data available

13. DISPOSAL CONSIDERATION

Waste material must be disposed of in accordance with the national and local regulations

14. TRANSPORT INFORMATION

DOT (US): Not regulated

UN proper shipping name:

Transport hazard class:

Packing group:

15. REGULATORY INFORMATION

US State Right to Know Laws:

California Proposition 65:

SPE Gel: This product can expose you to Barbiturates, which is known to the State of California to cause developmental harm.

For more information, go to: www.p65warnings.ca.gov.

International regulations:

EU Reach: No information available

16. OTHER INFORMATION

SDS Creation date: 8/28/2025

Revision: 1

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Helena Laboratories shall not be held liable for any damage resulting from handling or from contact with the above product.

HELENA LABORATORIES**SAFETY DATA SHEET**

DATE PREPARED: July 14, 2025

REVISION: 1

1. IDENTIFICATION**Product number:** 551151**Product identifier used on the label:** Acid Blue Stain**Recommended use of the chemical and restrictions in use:** For In-Vitro Diagnostic use.**Name, address, and telephone number of the chemical manufacturer, importer, or other responsible party:**

Helena Laboratories
1530 Lindbergh Dr.
Beaumont, TX 77707
Tel: (409) 842-3714
USA Toll Free 800-231-5663

Emergency phone number:

(409) 842-3714

2. HAZARD IDENTIFICATION**GHS Classification in accordance with the OSHA Hazard Communication Standard (29 CFR 1910.1200)**

Not a Hazardous Substance or Mixture

GHS Label elements

No signal word, no hazard statement(s), no precautionary statement(s) require

Other Hazards

None known

3. COMPOSITION/INFORMATION ON INGREDIENTS

Component Name	Chemical Name	CAS #	Concentration
Acid Blue Stain	Acid Blue 29	5850-35-1	100 %

This product contains no hazardous chemicals at or above Threshold limits

4. FIRST AID MEASURES

Description of first aid measures:

Inhalation: Remove person to fresh air. If not breathing, give artificial respiration. Consult a physician.

Skin contact: Immediately wash off with soap and plenty of water. If skin irritation occurs, get medical attention. Take off contaminated clothing.

Eye contact: Immediately flush eyes with plenty of water for at least 15 minutes. Consult a physician.

Ingestion: Do not induce vomiting. Rinse mouth with water. Consult a physician.

5. FIREFIGHTING MEASURES

Extinguishing media: Water spray, Dry chemical, carbon dioxide, or Alcohol resistant foam should be used.

Special hazards: None known

Special protective equipment and precautions for fire-fighters: Wear proper protective equipment and self contained breathing apparatus for firefighting if necessary.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment, and emergency procedures: Wear personal protective equipment. Soak up with inert absorbent material and transfer to container for proper disposal.

Methods and materials for containment and cleaning up: Do not allow material to enter drains/surface water/ground water. Ventilate the area.

7. HANDLING AND STORAGE

Precautions for safe handling: Wear personal protective equipment, avoid contact with skin and eye.

Conditions for safe storage, including any incompatibilities: Store as per specified instructions on label.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Exposure limits (AGGIH TLV & OSHA PEL):

Contains no substances with occupational exposure limit value.

Exposure controls

Respiratory protection: Not normally required

Skin protection: Wear appropriate protective gloves and suitable protective clothing.

Eye protection: Wear appropriate safety glasses or goggles.

Hygiene measures: Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Blue black powder

Odor: None

Odor threshold: No data available

pH: No data available

Melting point/Freezing point: No data available

Initial boiling point and range: No data available

Flash point: No data available

Evaporation rate: No data available

Flammability (solid, gas): No data available

Upper/lower flammability or explosive limits: No data available

Vapor pressure: No data available

Relative density: No data available

Solubility in water: Fully soluble

Partition coefficient: n-octanol/water: No data available

Auto-ignition temperature: No data available

Decomposition temperature: No data available

Viscosity: No data available

10. STABILITY AND REACTIVITY

Reactivity: Not normally reactive

Chemical stability: Stable under normal conditions

Possibility of hazard reaction: None known

Conditions to avoid: Heat, Hot surfaces, Flames

Incompatible materials: No data available

Hazard decomposition products: No data available

11. TOXICOLOGICAL INFORMATION**Potential health effects:**

Skin: May cause skin irritation.

Eye: May cause eye irritation

Inhalation: May be harmful if inhaled.

Ingestion: May be harmful if swallowed.

Carcinogenicity: No components are listed as carcinogen by AGGIH, IARC, OSHA or NTP.

Other important toxicological hazards: To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

12. ECOLOGICAL INFORMATION

Ecotoxicity: No data available

Persistence and degradability: No data available

Bioaccumulative potential: No data available

Mobility in soil: No data available

Other adverse effects: No data available

13. DISPOSAL CONSIDERATION

Observe all federal, state and local regulations.

14. TRANSPORT INFORMATION

DOT (US): Not regulated

UN number:

UN proper shipping name:

Transport hazard class:

Packing group:

15. REGULATORY INFORMATION

US State Right to Know Laws:

California Proposition 65: Not listed

EU Reach: No information available.

16. OTHER INFORMATION

SDS Creation date: July 14, 2025

Revision: 1

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Helena Laboratories shall not be held liable for any damage resulting from handling or from contact with the above product.

HELENA LABORATORIES**SAFETY DATA SHEET**

DATE PREPARED: 5/19/2025

REVISION: 1

1. IDENTIFICATION**Product identifier used on the label:** Citric Acid Destain**Product number:** 551959**Recommended use of the chemical and restrictions in use:** For In-Vitro Diagnostic use.**Name, address, and telephone number of the chemical manufacturer, importer, or other responsible party:**

Helena Laboratories
1530 Lindbergh Dr.
Beaumont, TX 77707
Tel: (409) 842-3714
USA Toll Free 800-231-5663

Emergency phone number:
(409) 842-3714**2. HAZARD IDENTIFICATION****Classification of the substrate or mixture****GHS Classification in accordance with OSHA Hazard Communication Standard (29 CFR 1910.1200)**

Eye Irritation Category 2A

Specific target organ toxicity – single exposure Category 3 (Respiratory system)

GHS Label elements**Signal word:** Warning**Hazard Statements:**

Cause serious eye irritation.

May cause respiratory irritation.

Precautionary Statements:

Avoid breathing dust.

Wear protective gloves / protective clothing / eye and face protection.

Wash skin thoroughly after handling.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Component Name	Chemical Name	CAS #	Concentration
Citric Acid Destain	Citric Acid	77-92-9	< = 100 %

4. FIRST AID MEASURES

Description of first aid measures:

Inhalation: Remove person to fresh air.

Skin contact: Take off contaminated clothing. Wash off with soap and plenty of water. If skin irritation occurs, get medical attention.

Eye contact: Immediately flush eyes with plenty of water. Remove contact lenses if present. Consult a physician.

Ingestion: If swallowed: immediately make victim drink water. Consult a physician.

5. FIREFIGHTING MEASURES

Suitable extinguishing media: Water spray, dry chemical, carbon dioxide , or foam.

Special hazards: Carbon Oxides. Combustible.

Special protective equipment and precautions for firefighters: Self-contained breathing apparatus (SCBA) and full protective clothing gear if necessary.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment, and emergency procedures: Wear personal protective equipment. Avoid inhalation of dust. Evacuate personnel to safe areas. Soak up with inert absorbent material and transfer to container for proper disposal.

Methods and materials for containment and cleaning up: Do not allow material to enter drains/surface water/ground water. Ventilate the area.

7. HANDLING AND STORAGE

Precautions for safe handling: Wear personal protective equipment, do not inhale, avoid contact with skin and eye.

Conditions for safe storage, including any incompatibilities: Keep in dry place.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Exposure limits (AGGIH TLV & OSHA PEL):

Contains no substances with occupational exposure limit values.

Exposure controls

Respiratory protection: No personal respiratory protective equipment normally required.
Use NIOSH approved respirators if necessary.

Skin protection: Wear appropriate protective gloves and suitable protective clothing.

Eye protection: Wear appropriate safety glasses, goggles or face shield.

Hygiene measures: Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Crystalline white powder

Odor: None

Odor Threshold: No data available

pH: No data available

Melting point/Freezing point: 153 – 159° C

Initial boiling point and range: Not applicable

Flash point: Not applicable

Evaporation Rate: No data available

Flammability (solid, gas): No data available

Upper/lower Flammability or explosive limits: Not applicable

Vapor pressure: < 0.1 hPa at 25° C

Relative density: No data available

Solubility in water: Soluble

Partition coefficient: n-octanol/water: No data available

Auto-ignition temperature: No data available

Decomposition temperature: No data available

Viscosity: No data available

10. STABILITY AND REACTIVITY

Reactivity: No decomposition if stored and applied as directed

Chemical stability: Stable under normal condition

Possibility of hazard reaction: Hazardous polymerization does not occur

Conditions to avoid: No information available

Incompatible materials: Oxidizing agents, strong Bases

Hazard decomposition products: In the event of fire: Carbon oxides

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure:

The substance or mixture has no acute oral toxicity

Carcinogenicity: Not listed as a carcinogen by IARC, OSHA or NTP.

Potential Chronic Health Effects: No known chronic health effects

Other important toxicological hazards: To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

12. ECOLOGICAL INFORMATION

Ecotoxicity: No data available

Persistence and degradability: Readily biodegradable.

Bioaccumulative potential: No data available

Mobility in soil: No data available.

Other adverse effects: No data available

13. DISPOSAL CONSIDERATION

Waste material must be disposed of in accordance with the national and local regulations

14. TRANSPORT INFORMATION

DOT (US): Not regulated

15. REGULATORY INFORMATION**US State Right to Know Laws:**

California Proposition 65: No Information available.

For more information, go to www.p65warnings.ca.gov

International regulations:

EU Reach: This substance has been registered according to Regulation (EC) No. 1907/2006 (REACH)

16. OTHER INFORMATION

SDS Creation date: 5/19/2025

Revision: 1

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Helena Laboratories shall not be held liable for any damage resulting from handling or from contact with the above product.